





FOREST 3

Forest 3
Rear wheel drive electronic wheelchair
with excellent outdoor capacities
and optimum indoor possibilities

















	Price	Image
Standard configuration		
S = Standard O = choice without extra cost □ = choice with extra cost - = not possible		

CONFIGURATION












Motor				
	6 km/h	S		 km/h
	10 km/h	□		
Seat width				
	400 mm	O		
	450 mm	O		
	500 mm	O		
Seat depth				
	420 mm	O		
	440 mm	O		
	460 mm	O		
	480 mm	O		
	500 mm	O		
	520 mm	O		
Seat height				
	415 mm	O		
	435 mm	O	ONLY IN COMBINATION WITH SE43 ELECTRICAL SEAT	
	465 mm	O		
	480 mm	O	ONLY IN COMBINATION WITH SE42 LIFT	
	485 mm	O	ONLY IN COMBINATION WITH SE43 ELECTRICAL SEAT	

Controller



R	Right	<input type="radio"/>		
L	Left	<input type="radio"/>		

Swing-away				
SE53	Swing-away controller	<input checked="" type="checkbox"/>		
SE53Z	Swing-away and fold down controller	<input type="checkbox"/>		
SE53A	Swing-away controller with 3 adjustable lock positions	<input type="checkbox"/>		
	Standard controller according to configuration	<input checked="" type="checkbox"/>		
LCD	Controller with integrated LCD touch screen	<input type="checkbox"/>		
SE64 - fixed	Remote actuator controls - fixed to armrest NOT POSSIBLE TO USE IN COMBINATION WITH SE53Z & SE53A	<input type="checkbox"/>		
SE64 - swivel	Remote actuator controls - on swivel of joystick NOT POSSIBLE TO USE IN COMBINATION WITH SE53Z & SE53A	<input type="checkbox"/>		
Controller options				
SE04	Chin control (equiped with LCD touch screen)	<input type="checkbox"/>		
SE08	Back control (controller positioned on backrest)	<input type="checkbox"/>		
SE09	Attendant control (2nd controller for attendant, positioned on backrest)	<input type="checkbox"/>		
SE11	Controller conv joystick ONLY IN COMBINATION WITH LE	<input type="radio"/>		
SE50	V-stab	<input type="checkbox"/>		
SE55	Reinforced joystick	<input type="checkbox"/>		
SE63	Softball joystick	<input type="checkbox"/>		
SE83	Mini joystick	<input type="checkbox"/>		
XD	Heavy duty power module (120A)	<input type="checkbox"/>		





ARMREST

	Armpad PU		S		
	Armpad Soft Touch comfort		<input type="checkbox"/>		
	Armpad artificial leather		<input type="checkbox"/>		
B66O	Gutted inclinable PU armrest NOT IN COMBINATION WITH SE74A MIRROR OR SE16	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66M	Hemiplegic arm rest with palm support NOT IN COMBINATION WITH B15 TABLE, SE74A MIRROR	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66ML	Hemiplegic arm rest with palm support - long model NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66MA	Hemiplegic arm rest with palm support - anatomic model NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66MAL	Hemiplegic arm rest with palm support - anatomic long model NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66OM	Hemiplegic arm rest with palm support - orientable NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66OML	Hemiplegic arm rest with palm support - long model - orientable NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66OMA	Hemiplegic arm rest with palm support - anatomic model - orientable NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		
B66OMAL	Hemiplegic arm rest with palm support - anatomic long model - orientable NOT IN COMBINATION WITH B15 TABLE, SE74A	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>		

FOOTREST

B06	Removable footrest		S		
BA6	Removable footrest short version (7cm shorter)		<input type="radio"/>		
BK6	Footrest with one footplate		<input type="checkbox"/>		




MANUAL LEGREST

BZ10	Monobloc footrest with legrest with side corners	<input type="checkbox"/>		
BZ7	Articulated length-correcting legrest, with articulated pallet (70mm longer)	<input type="checkbox"/>		
BZ8	Articulated length-correcting legrest, with articulated pallet	<input type="checkbox"/>		
BK7		<input type="checkbox"/>		
BL7	Amputee support	<input type="checkbox"/>		
BZE	Electrical articulated legrest	<input type="checkbox"/>		

COLORS (frame)

		Mat		Gloss			
C4	Bright red	○		○		-	
C5	Black	○		○		-	
C8	Yellow	○		○		-	
C29	Silver	○		○		-	
C53	Orange	○		○		-	
C55	Inox	○		○		-	
C56	Metallic blue	○		○		-	
C59	Champagne	○		○		-	
C71	Mirror blue	-		○		-	
C72	Mirror red	-		○		-	
C73	Mirror lime green	-		○		-	
C76	White	○		○		-	
C78	Mirror light blue	-		○		-	
C79	Soft blue	○		○		-	
C85	Bright green	○		○		-	
C87	Melon orange	○		○		-	
C88	Gold	○		○		-	
C89	Brown	○		○		-	
Additional colors with additional cost : contact your sales representative							

BACKREST

Type of backrest				
L34	Preshaped cushion			
	Soft Touch upholstery (grey)	<input type="radio"/>		
	3D Net upholstery (black)	<input type="checkbox"/>		
	Artificial leather upholstery (black)	<input type="checkbox"/>		
Tubular backrest : with preshaped cushion adjustable with alu clamps				
L34 Tub	Soft Touch upholstery (grey)	<input type="radio"/>		
	3D Net upholstery (black)	<input type="checkbox"/>		
	Artificial leather upholstery (black)	<input type="checkbox"/>		
	Tubular backrest with adjustable straps and black Nylon coating			
L44	Nylon upholstery	<input type="radio"/>		
	Soft Touch upholstery (grey)	<input type="radio"/>		
	3D Net upholstery (black)	<input type="checkbox"/>		
	Nylon PLUS upholstery	<input type="radio"/>		

SEAT

L35	Preshaped seat cushion			
	Soft Touch upholstery (grey)	<input type="radio"/>		
	3D Net upholstery (black)	<input type="checkbox"/>		
	3D Net + Nylon PLUS upholstery	<input type="checkbox"/>		
	Artificial leather upholstery (black)	<input type="checkbox"/>		
Preshaped seat cushion with velro				
L16F	No cushion (wooden plate)	<input type="radio"/>		
L16F+L35 Soft	Soft Touch upholstery (grey)	<input type="radio"/>		
L16F+L35 Net/Co	3D NET + Nylon PLUS	<input type="checkbox"/>		
L16F+L35 Skaï	Artificial leather upholstery (black)	<input type="checkbox"/>		
L16F+B23 nylon	Nylon upholstery	<input type="radio"/>		

SEAT & BACKREST ADJUSTMENT



SE30	Manual seat adjustment	<input type="radio"/>		
SE34	Manual backrest adjustment	<input type="radio"/>		
SE36	Manual backrest ; inclination to 30° with gas springs	<input type="checkbox"/>		
SE39	Electrical seat (20°) & electrical backrest (30°) inclination (2 actuators) (seat height w/o cushion: 415 or 465 mm)	<input type="checkbox"/>		
SE40	Electrical seat adjustment (20° - 1 actuator) (seat height w/o cushion: 415 or 465 mm)	<input type="checkbox"/>		
SE41	Electrical backrest adjustment (30° - 1 actuator) (seat height w/o cushion: 415 or 465 mm)	<input type="checkbox"/>		
SE42+SE40	LIFT of 300 mm including electrical seat adjustment (45°) (2 actuators) (seat height w/o cushion: 480 mm)	<input type="checkbox"/>		
SE42+SE39	LIFT of 300 mm including electrical seat (45°) & backrest (30°) adjustment (3 act.) (seat height w/o cushion: 480 mm)	<input type="checkbox"/>		
SE43	Electrical seat (50°) inclination (1 actuators) (seat height w/o cushion: 435 or 485 mm)	<input type="checkbox"/>		
SE43+SE41	Electrical seat (50°) & electrical backrest (30°) inclination (2 actuators) (seat height w/o cushion: 435 or 485 mm)	<input type="checkbox"/>		

WHEELS & TYRES



Front wheels

9"	(230 mm)			
	Pneumatic black tyre, black rim, line profile	<input type="radio"/>		
	Pneumatic grey tyre, black rim, line profile	<input type="radio"/>		
	Puncture-free black tyre, black rim, line profile	<input type="checkbox"/>		
10"	(260 mm)			
	Pneumatic black tyre, black rim, block profile	<input type="radio"/>		
	Pneumatic black tyre, black rim, line profile	<input type="radio"/>		
	Pneumatic grey tyre, black rim, block profile	<input type="checkbox"/>		
	Puncture-free black tyre, black rim, block profile	<input type="checkbox"/>		
	Puncture-free black tyre, black rim, line profile	<input type="checkbox"/>		
Rear wheels (3.00 x 8)				
	Pneumatic black tyre, black rim, block profile	<input type="radio"/>		
	Pneumatic black tyre, black rim, line profile	<input type="radio"/>		
	Pneumatic grey tyre, black rim, block profile	<input type="checkbox"/>		
	Puncture-free black tyre, black rim, block profile	<input type="checkbox"/>		














HEAD SUPORT SELECTION





L55	Headrest	<input type="checkbox"/>		
L58	Adjustable headrest in height and depth (stepless)	<input type="checkbox"/>		
L58C	Adjustable headrest in height, width, depth & rotation	<input type="checkbox"/>		






SEATING OPTIONS

SE16	Lateral side plates (in plexi) NOT IN COMBINATION WITH B66 / B660	<input type="checkbox"/>		
B22	Abductor	<input type="checkbox"/>		
L04	Side support	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>	
L04S	Side support (small)	<input type="radio"/> Right <input type="radio"/> Left	<input type="checkbox"/>	

GENERAL OPTIONS

B15	Plexi Table	<input type="checkbox"/>		
B15A	Adjustable table	<input type="checkbox"/>		
B31	Crutches holder NOT IN COMBINATION WITH SE42 & SE43	<input type="checkbox"/>		
B52	Serum holder	<input type="checkbox"/>		
B58	Safety belt	S		
SE01	Mudguards	<input type="checkbox"/>		
SE06	USB charger port	<input type="checkbox"/>		
SE18	Steering lock	<input type="checkbox"/>		
SE19	LED lights	<input type="checkbox"/>		
SE58	Taxi hooks	S		
SE24	Cupholder	<input type="checkbox"/>		
SE26	Pushing bar ONLY COMPATIBLE WITH L34 NOT COMPATIBLE WITH L34 Tub	<input type="checkbox"/>		
SE27	Top case NOT IN COMBINATION WITH TUBULAR BACKREST, SE42 & SE43 & L55 &	<input type="checkbox"/>		
SE28	Luggage rack NOT IN COMBINATION WITH SE42 & SE43	<input type="checkbox"/>		

SE37	Side support for femur	<input type="checkbox"/>		
SE74	Mirror (1 piece) NOT IN COMBINATION WITH B66 / B66O / B15 / B15A	<input type="checkbox"/>	<input type="radio"/> Right <input type="radio"/> Left	
SE74A	Swing-away mirror (1 piece) NOT IN COMBINATION WITH B66 / B66O / B15 / B15A WILL BE MOUNTED ON OPPOSITE SIDE OF	<input type="checkbox"/>	<input type="radio"/> Right <input type="radio"/> Left	
SE49	Bag Forest 3	<input type="checkbox"/>		

	Postural solutions - please refer to Neoflex to orderform
	Seating and positioning solutions - please refer to Physipro orderform
	Seating and positioning solutions - please refer to Vicair orderform
	Comfort gel applications - please refer to Gel cushions orderform
	Granulated Vacuum stabilisation solutions - please refer to Stabilo orderform

*Availability may differ depending on your region.
Please consult your sales manager for the most up-to-date information regarding all possibilities.*