


Forest 3/ Forest 3+ / Sigma / Timix

Turios




Material			
No.	Plan	Article	Quantity
1	Y024077	2080954	1
2	Charger	1916653	1
3*	Cable 1.5 m	1912299	1
4	M6 DIN985	1001680	2
5	M6x35 DIN7984	1903903	2
6	M8x10 ISO7380	1992177	1

* Not depicted

E001171 / A

1916727

SE06 LINX USB CHARGER				Plan No.
				E001171 / A
Scale: 1:2	Drafted	Verified	Confirmed	
Date	28/01/2021	28/01/2021	28/01/2021	
Initials	SD	SD	PJ	